

# Anti-CD19 (3G7) Mouse antibody



Catalog Number: 168195

## Product name

Anti-CD19 (3G7) Mouse antibody

## Specificity

Human

## Antibody description

CD19 (3G7) Mouse monoclonal antibody

## Formulation

Purified mouse monoclonal antibody in PBS(pH 7.4) containing with 0.03% Proclin300 and 50% glycerol.

## Storage

Store at 4°C short term. Store at -20°C long term.  
Avoid freeze / thaw cycle.

## Clonality

Monoclonal

## Ig Type

Mouse IgG1


## Applications

ICC/IF, FC


## Dilutions

ICC: 1/100

## Validations



Flow cytometric analysis of Jurkat T cells (black) and Ramos B cells (red), using anti-CD19 mouse mAb.



Confocal immunofluorescent analysis of Ramos (positive cell, left) and Jurkat (negative cell, right) using anti-CD19 mouse mAb (dilution 1:100).