

Anti-Ntrk2 antibody

Catalog Number: 176595

Product name

Anti-Ntrk2 antibody

Specificity

Mouse, Rat

Antibody description

Rabbit monoclonal antibody to Ntrk2

Preparation

This antigen of this antibody was recombinant protein within human trkb aa 100-450.

Formulation

Liquid, 1*TBS (pH7.4), 0.05% BSA, 40% Glycerol.
Preservative: 0.05% Sodium Azide.

Storage

Store at +4°C after thawing. Aliquot store at -20°C. Avoid repeated freeze / thaw cycles.

Clonality

Monoclonal

Ig Type

Rabbit IgG

Applications

WB

Dilutions

WB: 1:500-1:2,000

Validations




Fig1.; Western blot analysis of TrkB on different lysates. Proteins were transferred to a PVDF membrane and blocked with 5% BSA in PBS for 1 hour at room temperature. The primary antibody (1/500) was used in 5% BSA at room temperature for 2 hours. Goat Anti-Rabbit IgG - HRP Secondary Antibody (HA1001) at 1:5,000 dilution was used for 1 hour at room temperature.; Positive control; Lane 1: mouse brain tissue lysate; Lane 2: mouse cerebellum tissue lysate; Lane 3: rat cerebellum tissue lysate