

# USP19 antibody



Catalog Number: 116595

## Product name

USP19 antibody

## Immunogen

[Human USP19 Recombinant protein \(GST tag\)](#)

## Specificity

Human; other species not tested.

## Antibody description

USP19 Rabbit Polyclonal antibody. Positive IP detected in Jurkat cells. Positive WB detected in Jurkat cells. Observed molecular weight by Western-blot: 150 kDa

## Preparation

This antibody was obtained by immunization of USP19 recombinant protein (Accession Number: XM\_017005626). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IP


## Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IP: 1:200-1:2000

## Validations



Jurkat cells were subjected to SDS PAGE followed by western blot with Catalog No:116595(USP19 Antibody) at dilution of 1:600



IP Result of anti-USP19 (IP:Catalog No:116595, 4ug; Detection:Catalog No:116595 1:500) with Jurkat cells lysate 2400ug.