

Anti-P2Y12 Rabbit antibody



Catalog Number: 167595

Product name

Anti-P2Y12 Rabbit antibody

Specificity

Human, Mouse, Rat

Antibody description

P2Y12 Rabbit monoclonal antibody

Preparation

Antigen: A synthetic peptide of human P2Y12

Formulation

50nM Tris-Glycine(pH 7.4), 0.15M NaCl,
40%Glycerol, 0.01% sodium azide and 0.05% BSA

Storage

Store at 4°C short term. Store at -20°C long term.
Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type

Rabbit IgG

Applications


WB, IP

Dilutions

WB: 1/1000

IP: 1/20

Validations



Western blot detection of P2Y12 in A549,HL-60,U251,U87-MG cell lysates using P2Y12 Rabbit mAb(1:1000 diluted).Predicted band size:40kDa.Observed band size:52kDa.