

Anti-DYNLL1 Rabbit antibody



Catalog Number: 168595

Product name

Anti-DYNLL1 Rabbit antibody

Specificity

Human, Mouse, Rat

Antibody description

DYNLL1 Rabbit monoclonal antibody

Preparation

Antigen: A synthetic peptide of human DYNLL1

Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Storage

Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type

Rabbit IgG

Applications

WB, IHC-P, FC, IP

Dilutions


1:1000-1:5000

1:20

1:200

1:20

Validations



Western blot detection of DYNLL1 in 3T3,HeLa cell lysates using DYNLL1 Rabbit mAb(1:1000 diluted).Predicted band size:10kDa.Observed band size:10kDa.