

Human GPNMB/Osteoactivin (His Tag) recombinant protein



Catalog Number: 502995

General Information

Gene Name Synonym

Transmembrane glycoprotein HGFIN

Protein Construction

A DNA sequence encoding the human GPNMB isoform 2 (Q14956-2) extracellular domain (Met 1-Pro 474) was fused with a polyhistidine tag at the C-terminus.

Organism

Human

Expression Host

Human Cells

QC Testing

Purity

> 98 % as determined by SDS-PAGE

Endotoxin

< 1.0 EU per μg of the protein as determined by the LAL method

Stability

Samples are stable for up to twelve months from date of receipt at -70°C

Predicted N terminal

Ala 22

Molecular Mass

The recombinant human GPNMB consists of 464

amino acids and predicts a molecular mass of 52.1 kDa. In SDS-PAGE under reducing conditions, the apparent molecular mass of rhGPNMB is approximately 85-90 kDa due to glycosylation.

Formulation

Lyophilized from sterile PBS, pH 7.4

1. 5 % trehalose and mannitol are added as protectants before lyophilization.

2. Please contact us for any concerns or special requirements.

Usage Guide

Storage


Store it under sterile conditions at -20°C to -80°C .

It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

Reconstitution

Adding sterile water, prepare a stock solution of 0.25 mg/ml. Concentration is measured by UV-Vis.

SDS-PAGE



Human GPNMB / Osteoactivin Protein (His Tag)
SDS-PAGE