

POGK antibody



Catalog Number: 113995

Product name

POGK antibody

Specificity

Human; other species not tested.

Antibody description

POGK Rabbit Polyclonal antibody. Positive IF detected in HeLa cells. Positive WB detected in A549 cells, human testis tissue, MCF-7 cells. Observed molecular weight by Western-blot: 85 kDa

Preparation

This antibody was obtained by immunization of POGK recombinant protein (Accession Number: XM_005245370). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, IF, WB


Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IF: 1:10-1:100

Validations



A549 cells were subjected to SDS PAGE followed by western blot with Catalog No:113995(POGK Antibody) at dilution of 1:600



Immunofluorescent analysis of HeLa cells using Catalog No:113995(POGK Antibody) at dilution of 1:25 and Rhodamine-Goat anti-Rabbit IgG