

SCGN antibody

Catalog Number: 114995

Product name

SCGN antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

SCGN Rabbit Polyclonal antibody. Positive IP detected in mouse brain tissue. Positive WB detected in mouse liver tissue, fetal human brain tissue, mouse colon tissue, mouse pancreas tissue, mouse stomach tissue, rat brain tissue. Observed molecular weight by Western-blot: 32 kDa

Preparation

This antibody was obtained by immunization of SCGN recombinant protein (Accession Number: NM_006998). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IP


Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IP: 1:200-1:2000

Validations



mouse liver tissue were subjected to SDS PAGE followed by western blot with Catalog No:114995(SCGN antibody) at dilution of 1:300



IP Result of anti-SCGN (IP:Catalog No:114995, 4ug; Detection:Catalog No:114995 1:500) with mouse brain tissue lysate 4000ug.