

# Anti-Fyn Rabbit antibody

Catalog Number: 166995

## Product name

Anti-Fyn Rabbit antibody

## Specificity

Human, Mouse, Rat

## Antibody description

Fyn Rabbit polyclonal antibody

## Preparation

Antigen: CTATEPQYQPGENL

## Formulation

10mM PBS with 0.02% sodium azide 50% glycerol.

## Storage

Store at 4°C short term. Store at -20°C long term.  
Avoid freeze / thaw cycle.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG


## Applications

WB

## Dilutions

WB: 1/1000

## Validations



Western blot analysis of extracts from Ramos, Jurkat, C6, 3T3 and HeLa cell lysates using Fyn Rabbit pAb (1:1000 diluted). Predicted band size: 59kDa. Observed band size: 45-59kDa.