

Anti-Bcl G/BCL2L14 Rabbit antibody



Catalog Number: 167995

Product name

Anti-Bcl G/BCL2L14 Rabbit antibody

Specificity

Human, Mouse, Rat

Antibody description

Bcl G/BCL2L14 Rabbit monoclonal antibody

Preparation

Antigen: Recombinant protein of human Bcl G/BCL2L14

Formulation

50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA

Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type

Rabbit IgG

Applications

WB, IHC-P, ICC/IF, FC, IP

Dilutions

WB: 1/1000


IHC: 1/20

ICC/IF: 1/50

FC: 1/20

IP: 1/20

Validations



Western blot detection of Bcl G/BCL2L14 in HeLa,A549,HL-60,U251,U87-MG cell lysates using Bcl G/BCL2L14 Rabbit mAb(1:1000 diluted).Predicted band size:37kDa.Observed band size:43kDa.