

Anti-CD32B Rabbit antibody



Catalog Number: 168995

Product name

Anti-CD32B Rabbit antibody

Specificity

Human

Antibody description

CD32B Rabbit monoclonal antibody

Preparation

Antigen: A synthetic peptide of human CD32B

Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Storage

Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type

Rabbit IgG

Applications


WB, IHC-P

Dilutions

1:2000-1:10000

1:200

Validations



Western blot detection of CD32B in C6 cell lysates using CD32B Rabbit mAb(1:1000 diluted).Predicted band size:34kDa.Observed band size:44kDa.