

Mouse PILRB (Fc Tag) recombinant protein



Catalog Number: 501295

General Information

Gene Name Synonym

Activating receptor PILR-beta; Cell surface receptor FDFACT

Protein Construction

A DNA sequence encoding the mouse PILRB1 (Q2YFS2) extracellular domain (Met 1-Gly 193) was fused with the Fc region of human IgG1 at the C-terminus.

Organism

Mouse

Expression Host

Human Cells

QC Testing

Purity

> 85 % as determined by SDS-PAGE

Endotoxin

< 1.0 EU per μg of the protein as determined by the LAL method

Stability

Samples are stable for up to twelve months from date of receipt at -70°C

Predicted N terminal

Gly 29

Molecular Mass

The secreted recombinant mouse PILRB1/Fc is a

disulfide-linked homodimer. The reduced monomer comprises 406 amino acids and has a calculated molecular mass of 45.8 kDa. As a result of glycosylation, the apparent molecular mass of rmIL1RAPL1/Fc monomer is approximately 55-60 kDa in SDS-PAGE under reducing conditions.

Formulation

Lyophilized from sterile PBS, pH 7.4

1. 5 % trehalose and mannitol are added as protectants before lyophilization.

2. Please contact us for any concerns or special requirements.

Usage Guide

Storage


Store it under sterile conditions at -20°C to -80°C .

It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

Reconstitution

Adding sterile water, prepare a stock solution of 0.25 mg/ml. Concentration is measured by UV-Vis.

SDS-PAGE



Mouse PILRB1 Protein (Fc Tag) SDS-PAGE