

# NFYB Polyclonal Antibody



Catalog Number: 165295

## Product name

NFYB Polyclonal Antibody

## Specificity

Human, Mouse

## Antibody description

Polyclonal antibody to NFYB

## Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 1-207 of human NFYB (NP\_006157.1).

## Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Storage

Store at -20°C. Avoid freeze / thaw cycles.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


WB, IF

## Dilutions

WB 1:500 - 1:2000


IF 1:50 - 1:100

## Validations



Western blot - NFYB Polyclonal Antibody

Western blot analysis of extracts of various cell lines, using NFYB antibody at 1:500 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL



Basic Kit .Exposure time: 90s.

Immunofluorescence - NFYB Polyclonal Antibody

Immunofluorescence analysis of A-549 cells using NFYB antibody . Blue: DAPI for nuclear staining.