

# FOXD4L6 antibody



Catalog Number: 110695

---

## Product name

FOXD4L6 antibody

## Specificity

Human; other species not tested.

## Antibody description

FOXD4L6 Rabbit Polyclonal antibody. Positive WB detected in A431 cells, HepG2 cells, Jurkat cells. Observed molecular weight by Western-blot: 46 kDa

## Preparation

This antibody was obtained by immunization of FOXD4L6 recombinant protein (Accession Number: NM\_001085476). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:500-1:5000

## Validations



A431 cells were subjected to SDS PAGE followed by western blot with Catalog No:110695(FOXD4L6 antibody) at dilution of 1:1000