

Anti-Angptl4 antibody



Catalog Number: 176695

Product name

Anti-Angptl4 antibody

Specificity

Mouse

Antibody description

Rabbit polyclonal antibody to Angptl4

Preparation

This antigen of this antibody was synthetic peptide within human angiopoietin-like 4 aa 60-122.

Formulation

Liquid, 1*PBS (pH7.4), 0.2% BSA, 40% Glycerol.
Preservative: 0.05% Sodium Azide.

Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB

Dilutions

WB: 1:500

Validations




Fig1: Western blot analysis on mouse liver tissue lysates using anti-angiopoietin-like 4 rabbit polyclonal antibodies.