

Anti-Desmocollin 3 Rabbit antibody



Catalog Number: 168695

Product name

Rabbit IgG

Anti-Desmocollin 3 Rabbit antibody

Applications

Specificity

WB, ICC/IF, FC

Human, Mouse, Rat

Dilutions

Antibody description

1:2000-1:10000

Desmocollin 3 Rabbit monoclonal antibody

1:20-1:50

Preparation


1:20

Antigen: A synthetic peptide of human
Desmocollin 3

Validations

Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M
NaCl, 40%Glycerol, 0.01% sodium azide and
0.05% BSA.



Storage

Store at 4°C short term. Aliquot and store at -20°C
long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal

Western blot detection of Desmocollin 3 in
K562, HeLa cell lysates using Desmocollin 3 Rabbit
mAb(1:1000 diluted). Predicted band
size:100kDa. Observed band size:100kDa.

Ig Type