

HEBP2 antibody

Catalog Number: 111395

Product name

HEBP2 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

HEBP2 Rabbit Polyclonal antibody. Positive WB detected in human brain tissue, human kidney tissue. Positive IP detected in HEK-293 cells. Observed molecular weight by Western-blot: 23kd

Preparation

This antibody was obtained by immunization of HEBP2 recombinant protein (Accession Number: NM_014320). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IP


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IP: 1:200-1:2000

Validations



human brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:111395(HEBP2 antibody) at dilution of 1:500



IP Result of anti-HEBP2 (IP:Catalog No:111395, 3ug; Detection:Catalog No:111395 1:500) with HEK-293 cells lysate 2400ug.