

GLA Polyclonal Antibody



Catalog Number: 163395

Product name

GLA Polyclonal Antibody

Specificity

Human, Mouse

Antibody description

Polyclonal antibody to GLA

Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 150-429 of human GLA (NP_000160.1).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB, IHC, IF


Dilutions

WB 1:500 - 1:2000

IHC 1:50 - 1:200


IF 1:50 - 1:100

Validations




Western blot - GLA Polyclonal Antibody

Western blot analysis of extracts of various cell lines, using GLA antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure



Immunohistochemistry - GLA Polyclonal Antibody

Immunohistochemistry of paraffin-embedded human liver using GLA antibody at dilution of




GLA Polyclonal Antibody

Catalog Number: 163395

Immunohistochemistry - GLA Polyclonal Antibody

Immunohistochemistry of paraffin-embedded mouse kidney using GLA antibody at dilution of



1:100 (40x lens).

Immunofluorescence - GLA Polyclonal Antibody

Immunofluorescence analysis of HeLa cells using GLA antibody . Blue: DAPI for nuclear staining.