

NT5DC1 antibody



Catalog Number: 113395

Product name

NT5DC1 antibody

Specificity

Human; other species not tested.

Antibody description

NT5DC1 Rabbit Polyclonal antibody. Positive WB detected in MCF7 cells, HepG2 cells, Jurkat cells. Positive IF detected in HepG2 cells. Observed molecular weight by Western-blot: 52 kDa

Preparation

This antibody was obtained by immunization of NT5DC1 recombinant protein (Accession Number: BC015138). Purification method: Antigen Affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IF


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IF: 1:50-1:500

Validations



MCF7 cells were subjected to SDS PAGE followed by western blot with Catalog No:113395(NT5DC1 antibody) at dilution of 1:1000



Immunofluorescent analysis of (10% Formaldehyde) fixed HepG2 cells using Catalog No:113395(NT5DC1 Antibody) at dilution of 1:50 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)