

CLMP antibody



Catalog Number: 109395

Product name

CLMP antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

CLMP Rabbit Polyclonal antibody. Positive WB detected in mouse brain tissue, 3T3-L1 cells, rat brain tissue. Observed molecular weight by Western-blot: 45 kDa

Preparation

This antibody was obtained by immunization of CLMP recombinant protein (Accession Number: NM_024769). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:2000

Validations



mouse brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:109395(ASAM Antibody) at dilution of 1:600