

# Anti-alpha 1 Antitrypsin Rabbit antibody



Catalog Number: 167795

## Product name

Anti-alpha 1 Antitrypsin Rabbit antibody

## Specificity

Human, Mouse, Rat

## Antibody description

alpha 1 Antitrypsin Rabbit monoclonal antibody

## Preparation

Antigen: A synthetic peptide of human alpha 1 Antitrypsin

## Formulation

50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA

## Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

## Clonality

Monoclonal

## Ig Type

Rabbit IgG

## Applications


WB, IP

## Dilutions

WB: 1/1000

IP: 1/20

## Validations



Western blot detection of alpha 1 Antitrypsin in HeLa, A549, MCF-7, U251, HL-60 cell lysates using alpha 1 Antitrypsin Rabbit mAb(1:1000 diluted). Predicted band size: 46kDa. Observed band size: 55kDa.