

ADM2 antibody

Catalog Number: 107795

Product name

ADM2 antibody

Specificity

Human, Mouse; other species not tested.

Antibody description

ADM2 Rabbit Polyclonal antibody. Positive IHC detected in human stomach tissue, human heart tissue, mouse heart tissue, mouse stomach tissue.

Preparation

This antibody was obtained by immunization of ADM2 recombinant protein (Accession Number: BC012864). Purification method: Antigen Affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, IHC

Dilutions


Recommended Dilution:

IHC: 1:20-1:200

Validations



Immunohistochemistry of paraffin-embedded human heart tissue slide using Catalog No:107795(ADM2 Antibody) at dilution of 1:50 under 40x lens



Immunohistochemistry of paraffin-embedded human stomach tissue slide using Catalog No:107795(ADM2 Antibody) at dilution of 1:50 under 40x lens