

Anti-HSPA2 Rabbit antibody



Catalog Number: 168795

Product name

Anti-HSPA2 Rabbit antibody

Specificity

Human, Mouse, Rat

Antibody description

HSPA2 Rabbit monoclonal antibody

Preparation

Antigen: Recombinant protein of human HSPA2

Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Storage

Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type

Rabbit IgG

Applications


WB, IHC-P

Dilutions

1:1000-1:5000

1:20-1:100

Validations



Western blot detection of HSPA2 in Rat Brain, C6, 3T3 cell lysates using HSPA2 Rabbit mAb(1:1000 diluted). Predicted band size: 70kDa. Observed band size: 70kDa.