

Catalog Number: 100006

Product name

Anti-TNFR1/CD120a/TNFRSF1A antibody

Immunogen

[Human TNFR1/CD120a/TNFRSF1A \(His Tag\) recombinant protein](#)

Specificity

Human TNFR1 / CD120a / TNFRSF1A

Antibody description

Rabbit monoclonal to TNFR1/CD120a/TNFRSF1A

Preparation

This antibody was obtained from a rabbit immunized with purified, recombinant Human TNFR1 / CD120a / TNFRSF1A (rh TNFR1 / CD120a / TNFRSF1A; NP_001056.1; Met1-Thr211).

Formulation

0.2 µm filtered solution in PBS

Storage

This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Preservative-Free.

Sodium azide is recommended to avoid contamination (final concentration 0.05%-0.1%). It is toxic to cells and should be disposed of properly. Avoid repeated freeze-thaw cycles.

Clonality

Monoclonal

Ig Type

Rabbit IgG


Applications

IF, ICC/IF

Dilutions

ICC/IF: 10-25 µg/mL

Validations



TNFR1 / CD120a / TNFRSF1A Antibody, Rabbit MAb, Immunofluorescence

Immunofluorescence staining of Human TNFR1 in Hela cells. Cells were fixed with 4% PFA, permeabilized with 0.3% Triton X-100 in PBS, blocked with 10% serum, and incubated with rabbit anti-Human TNFR1 monoclonal antibody (15 µg/ml) at 4°C overnight. Then cells were stained with the Alexa Fluor® 488-conjugated Goat Anti-rabbit IgG secondary antibody (green) and counterstained with DAPI (blue). Positive staining was localized to cytoplasm.