

DMRTA1 antibody



Catalog Number: 110006

Product name

DMRTA1 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

DMRTA1 Rabbit Polyclonal antibody. Positive WB detected in mouse liver tissue. Positive IP detected in mouse liver tissue. Observed molecular weight by Western-blot: 60-70 kDa

Preparation

This antibody was obtained by immunization of DMRTA1 recombinant protein (Accession Number: NM_022160). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, IP, WB


Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IP: 1:200-1:2000

Validations



mouse liver tissue were subjected to SDS PAGE followed by western blot with Catalog No:110006(DMRTA1 Antibody) at dilution of 1:600



IP Result of anti-DMRTA1 (IP:Catalog No:110006, 3ug; Detection:Catalog No:110006 1:500) with mouse liver tissue lysate 3000ug.