

# GNMT antibody

Catalog Number: 111006

## Product name

GNMT antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

GNMT Rabbit Polyclonal antibody. Positive WB detected in mouse liver tissue, mouse pancreas tissue. Positive IP detected in mouse liver tissue. Observed molecular weight by Western-blot: 33 kDa

## Preparation

This antibody was obtained by immunization of GNMT recombinant protein (Accession Number: NM\_018960). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IP


## Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IP: 1:200-1:2000


## Validations



mouse liver tissue were subjected to SDS PAGE followed by western blot with Catalog No:111006(GNMT antibody) at dilution of 1:500



mouse liver tissue were subjected to SDS PAGE followed by western blot with Catalog No:111006(GNMT Antibody) at dilution of 1:600



IP Result of anti-GNMT (IP:Catalog No:111006, 4ug; Detection:Catalog No:111006 1:600) with mouse liver tissue lysate 6800ug.