

PDK4 Polyclonal Antibody



Catalog Number: 164006

Product name

Rabbit IgG

PDK4 Polyclonal Antibody

Applications**Specificity**

WB, IF

Human, Mouse

Dilutions**Antibody description**


WB 1:500 - 1:2000

Polyclonal antibody to PDK4

IF 1:10 - 1:100

Preparation**Validations**

Antigen: A synthetic peptide of human PDK4

**Formulation**

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Immunofluorescence - PDK4 Polyclonal Antibody

Polyclonal

Immunofluorescence analysis of U2OS cells using PDK4 antibody . Blue: DAPI for nuclear staining.

Ig Type