

ERC2 antibody



Catalog Number: 110406

Product name

ERC2 antibody

Immunogen

[Human ERC2 Recombinant protein \(GST tag\)](#)

Specificity

human, mouse, rat; other species not tested.

Antibody description

ERC2 Rabbit Polyclonal antibody. Positive IF detected in Hela. Positive WB detected in human brain tissue, rat brain tissue. Observed molecular weight by Western-blot: 115-120kd

Preparation

This antibody was obtained by immunization of ERC2 recombinant protein (Accession Number: XM_017006158). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IF


Dilutions

Recommended Dilution:


WB: 1:1000-1:10000

IF: 1:10-1:100

Validations



human brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:110406(ERC2 antibody) at dilution of 1:500



Immunofluorescent analysis of Hela cells, using ERC2 antibody Catalog No:110406 at 1:25 dilution and Rhodamine-labeled goat anti-rabbit IgG (red).