

Anti-SLC6A17 antibody



Catalog Number: 176406

Product name

Anti-SLC6A17 antibody

Ig Type

Rabbit IgG

Specificity

Human, Mouse, Rat

Applications

WB

Antibody description

Rabbit polyclonal antibody to SLC6A17

Dilutions

WB: 1:500-1:2000

Preparation

This antigen of this antibody was synthesized peptide derived from ntt4 . at aa range: 240-320

Formulation

Liquid, Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

Clonality

Polyclonal

Validations




Fig1:; Western Blot analysis of extracts from rat stomach, using NTT4 Polyclonal Antibody.. Secondary antibody (catalog#: HA1001) was diluted at 1:20000