

# GABPB2 antibody



Catalog Number: 110806

---

## Product name

GABPB2 antibody

## Specificity

human, mouse, rat; other species not tested.

## Antibody description

GABPB2 Rabbit Polyclonal antibody. Positive WB detected in Jurkat cells, HEK-293 cells, HeLa cells. Observed molecular weight by Western-blot: 49kd

## Preparation

This antibody was obtained by immunization of GABPB2 recombinant protein (Accession Number: XM\_017000249). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:200-1:2000

## Validations



Jurkat cells were subjected to SDS PAGE followed by western blot with Catalog No:110806(GABPB2 antibody) at dilution of 1:500