

# Human Tie2/CD202b/TEK (His Tag) recombinant protein



Catalog Number: 503806

## General Information

### Gene Name Synonym

Endothelial tyrosine kinase; Tunica interna endothelial cell kinase; Tyrosine kinase with Ig and EGF homology domains-2; Tyrosine-protein kinase receptor TEK; Tyrosine-protein kinase receptor TIE-2; p140 TEK

### Protein Construction

A DNA sequence encoding the extracellular domain (Met 1-Lys 745) of human Tie2 (NP\_000450.2) precursor was fused with a polyhistidine tag at the C-terminus.

### Organism

Human

### Expression Host

Human Cells

## QC Testing

### Activity

Measured by its binding ability in a functional ELISA. Immobilized recombinant human Tie2 at 2 µg/ml (100 µl/well) can bind human Angiopoietin-2 at a linear range of 1.28-160 ng/ml.

### Purity

> 95 % as determined by SDS-PAGE

### Endotoxin

< 1.0 EU per µg of the protein as determined by the LAL method

### Stability

Samples are stable for up to twelve months from date of receipt at -70°C

### Predicted N terminal

Ala 23

### Molecular Mass

The recombinant human Tie2 comprises 734 amino acids and predicts a molecular mass of 82 kDa. As a result of glycosylation, the apparent molecular mass of rhTie2 is approximately 95-105 kDa in SDS-PAGE under reducing conditions.

### Formulation

Lyophilized from sterile PBS, pH 7.4

1. 5 % trehalose and mannitol are added as protectants before lyophilization.

2. Please contact us for any concerns or special requirements.

## Usage Guide


### Storage

Store it under sterile conditions at -20°C to -80°C. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

### Reconstitution

Adding sterile water, prepare a stock solution of 0.25 mg/ml. Concentration is measured by UV-Vis.

### SDS-PAGE



Human Tie2 / CD202b / TEK Protein (His Tag)  
SDS-PAGE