

PLEK Polyclonal Antibody



Catalog Number: 164806

Product name

PLEK Polyclonal Antibody

Specificity

Human, Mouse, Rat

Antibody description

Polyclonal antibody to PLEK

Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 1-350 of human PLEK (NP_002655.2).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


WB, IHC

Dilutions

WB 1:500 - 1:2000


IHC 1:50 - 1:200

Validations




Western blot - PLEK Polyclonal Antibody

Western blot analysis of extracts of various cell lines, using PLEK antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 10s.



Immunohistochemistry - PLEK Polyclonal Antibody

Immunohistochemistry of paraffin-embedded rat heart using PLEK antibody at dilution of 1:200



(40x lens).

PLEK Polyclonal Antibody



Catalog Number: 164806

Immunohistochemistry - PLEK Polyclonal Antibody human stomach using PLEK antibody at dilution of 1:200 (40x lens).
Immunohistochemistry of paraffin-embedded