

# ABCF2 antibody



Catalog Number: 107806

## Product name

ABCF2 antibody

## Immunogen

[Human ABCF2 Recombinant protein \(GST tag\)](#)

## Specificity

Human, Mouse; other species not tested.

## Antibody description

ABCF2 Rabbit Polyclonal antibody. Positive WB detected in HEK-293 cells, HeLa cells, NIH/3T3 cells. Positive IHC detected in human colon cancer tissue. Observed molecular weight by Western-blot: 71 kDa

## Preparation

This antibody was obtained by immunization of ABCF2 recombinant protein (Accession Number: XM\_017011640). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IHC


## Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IHC: 1:20-1:200

## Validations



HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:107806(ABCF2 antibody) at dilution of 1:1000



Immunohistochemical of paraffin-embedded human colon cancer using Catalog No:107806(ABCF2 antibody) at dilution of 1:100 (under 25x lens)