

# C9orf9 antibody



Catalog Number: 108806

## Product name

C9orf9 antibody

## Immunogen

[Human C9orf9 Recombinant protein \(GST tag\)](#)

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

C9orf9 Rabbit Polyclonal antibody. Positive WB detected in mouse testis tissue, mouse liver tissue, rat testis tissue. Positive IF detected in A431 cells. Observed molecular weight by Western-blot: 18, 30 kDa

## Preparation

This antibody was obtained by immunization of C9orf9 recombinant protein (Accession Number: NM\_001316897). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IF


## Dilutions

Recommended Dilution:


WB: 1:200-1:1000

IF: 1:10-1:100

## Validations



mouse testis tissue were subjected to SDS PAGE followed by western blot with Catalog No:108806(C9orf9 Antibody) at dilution of 1:300



Immunofluorescent analysis of A431 cells using Catalog No:108806(C9orf9 Antibody) at dilution of 1:25 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)