

DNAJC16 antibody



Catalog Number: 110106

Product name

DNAJC16 antibody

Immunogen

[Human DNAJC16 Recombinant protein \(His tag\)](#)

Specificity

human, mouse, rat; other species not tested.

Antibody description

DNAJC16 Rabbit Polyclonal antibody. Positive WB detected in mouse brain tissue. Observed molecular weight by Western-blot: 91kd

Preparation

This antibody was obtained by immunization of DNAJC16 recombinant protein (Accession Number: NM_015291). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

ELISA, WB

Dilutions

Recommended Dilution:
WB: 1:500-1:5000

Validations



mouse brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:110106(DNAJC16 antibody) at dilution of 1:1500