

PPP2R5D antibody



Catalog Number: 114106

Product name

PPP2R5D antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

PPP2R5D Rabbit Polyclonal antibody. Positive IP detected in Y79 cells. Positive WB detected in Y79 cells, human brain tissue. Observed molecular weight by Western-blot: 70 kDa

Preparation

This antibody was obtained by immunization of PPP2R5D recombinant protein (Accession Number: NM_006245). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IP


Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IP: 1:200-1:1000

Validations



Y79 cells were subjected to SDS PAGE followed by western blot with Catalog No:114106(PPP2R5D antibody) at dilution of 1:400



IP Result of anti-PPP2R5D (IP:Catalog No:114106, 4ug; Detection:Catalog No:114106 1:400) with Y79 cells lysate 2000ug.