

Anti-Scn10a antibody



Catalog Number: 175106

Product name

Anti-Scn10a antibody

Specificity

Mouse, Rat

Antibody description

Rabbit polyclonal antibody to Scn10a

Preparation

This antigen of this antibody was recombinant protein

Formulation

Liquid, 1*PBS (pH7.4), 0.2% BSA, 50% Glycerol.
Preservative: 0.05% Sodium Azide.

Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB, IHC-P

Dilutions

WB: 1:100-1:500

IHC-P: 1:50-1:200

Validations





Fig1: Western blot analysis of NaV1.8 on rat spinal cord lysate using anti-NaV1.8 antibody at



1/100 dilution.

Fig2: Immunohistochemical analysis of paraffin-embedded mouse brain tissue using anti-NaV1.8 antibody. Counter stained with hematoxylin.