

Anti-DsbA (1A2) Mouse antibody



Catalog Number: 169106

Product name

Anti-DsbA (1A2) Mouse antibody

Specificity

E.coli

Antibody description

DsbA (1A2) Mouse monoclonal antibody

Preparation

Antigen: Purified recombinant E.coli. DsbA protein fragments over expressed in E.coli.

Formulation

Purified mouse monoclonal in buffer containing 0.1M Tris-Glycine(pH 7.4 150 mM NaCl)with 0.02% sodium azide 0.1mg/mlBSA and 50% glycerol.

Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type

Mouse IgG2b


Applications

WB


Dilutions

WB: 1/1000

Validations



Western blot detection of DSBA in AKT1,AKT2,AKT3 and DSBA recombinant antigen fragments using DSBA mouse mAb (1:1000



diluted). 19KDa. Observed band size: 24KDa. Predicted band size: 24KDa.