

# RAC2 antibody

Catalog Number: 114506

## Product name

RAC2 antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

RAC2 Rabbit Polyclonal antibody. Positive WB detected in HL-60 cells. Positive IP detected in HL-60 cells. Observed molecular weight by Western-blot: 21 kDa

## Preparation

This antibody was obtained by immunization of RAC2 recombinant protein (Accession Number: NM\_002872). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IP


## Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IP: 1:200-1:1000

## Validations



HL-60 cells were subjected to SDS PAGE followed by western blot with Catalog No:114506(RAC2 Antibody) at dilution of 1:600



IP Result of anti-RAC2 (IP:Catalog No:114506, 3ug; Detection:Catalog No:114506 1:300) with HL-60 cells lysate 4000ug.