

RISC antibody

Catalog Number: 107506

Product name

RISC antibody

Specificity

Human, Mouse; other species not tested.

Antibody description

RISC Mouse Monoclonal antibody. Positive WB detected in human testis tissue, human heart tissue, MCF7 cells. Positive IF detected in C6 cells. Observed molecular weight by Western-blot: 35-40kd/50-60kd

Preparation

This antibody was obtained by immunization of RISC recombinant protein (Accession Number: AF330052). Purification method: Caprylic Acid/Ammonium Sulfate Precipitation.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Monoclonal

Ig Type

Mouse IgG2b

Applications

ELISA, WB, IF


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IF: 1:10-1:100

Validations



human testis tissue were subjected to SDS PAGE followed by western blot with Catalog No:107506(Scepep1 antibody) at dilution of 1:1000



Immunofluorescent analysis of C6 cells using Catalog No:107506(Scepep1 Antibody) at dilution of 1:25 and Rhodamine-Goat anti-Mouse IgG