

Anti-HSP40 Mouse antibody



Catalog Number: 168506

Product name

Anti-HSP40 Mouse antibody

Specificity

Human

Antibody description

HSP40 Mouse monoclonal antibody

Preparation

Antigen: Recombinant Protein of Hsp40

Formulation

PBS, pH 7.4, containing 0.02% sodium azide as Preservative and 50% Glycerol.

Storage

Store at -20°C

Clonality

Monoclonal

Ig Type

Mouse IgG


Applications

WB

Dilutions

1:1000

Validations



Western blot detection of HSP40 in Jurkat cell lysates using HSP40 Mouse mAb(1:1000 diluted). Predicted band size:40KDa. Observed band size:40KDa.