

GCLC antibody

Catalog Number: 110906

Product name

GCLC antibody

IP: 1:200-1:2000

IHC: 1:20-1:200


Specificity

Human, Mouse, Rat; other species not tested.

Validations

Antibody description

GCLC Rabbit Polyclonal antibody. Positive WB detected in Jurkat cells, HEK-293 cells, L02 cells, mouse kidney tissue, mouse lung tissue. Positive IP detected in mouse kidney tissue. Positive IHC detected in human kidney tissue. Observed molecular weight by Western-blot: 73kd




Preparation

This antibody was obtained by immunization of GCLC recombinant protein (Accession Number: BC022487). Purification method: Antigen affinity purified.

Jurkat cells were subjected to SDS PAGE followed by western blot with Catalog No:110906(GCLC antibody) at dilution of 1:200

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.



Storage

Store at -20°C. DO NOT ALIQUOT


Clonality

Polyclonal

Immunohistochemical of paraffin-embedded human kidney using Catalog No:110906(GCLC antibody) at dilution of 1:50 (under 10x lens)

Ig Type

Rabbit IgG



Applications

ELISA, WB, IHC, IP

Dilutions

Recommended Dilution:


WB: 1:200-1:2000

Immunohistochemical of paraffin-embedded human kidney using Catalog No:110906(GCLC antibody) at dilution of 1:50 (under 10x lens)

GCLC antibody

Catalog Number: 110906

antibody) at dilution of 1:50 (under 40x lens)



IP Result of anti-GCLC (IP:Catalog No:110906, 3ug; Detection:Catalog No:110906 1:500) with mouse kidney tissue lysate 4000ug.