

TCF7 antibody



Catalog Number: 115906

Product name

TCF7 antibody

Specificity

Human, Mouse; other species not tested.

Antibody description

TCF7 Rabbit Polyclonal antibody. Positive WB detected in mouse thymus tissue. Observed molecular weight by Western-blot: 43 kDa

Preparation

This antibody was obtained by immunization of TCF7 recombinant protein (Accession Number: NM_003202). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:2000

Validations



mouse thymus tissue were subjected to SDS PAGE followed by western blot with Catalog No:115906(TCF7 Antibody) at dilution of 1:600