

ZBP1 antibody

Catalog Number: 116906

Product name

ZBP1 antibody

Specificity

human; other species not tested.

Antibody description

ZBP1 Rabbit Polyclonal antibody. Positive IHC detected in human liver cancer. Positive WB detected in Jurkat cells, HL-60 cells. Observed molecular weight by Western-blot: 68kd

Preparation

This antibody was obtained by immunization of ZBP1 recombinant protein (Accession Number: BC028218). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IHC: 1:20-1:200


Validations



Jurkat cells were subjected to SDS PAGE followed by western blot with Catalog No:116906(ZBP1 antibody) at dilution of 1:400



Immunohistochemical of paraffin-embedded human liver cancer using Catalog No:116906(ZBP1 antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemical of paraffin-embedded human liver cancer using Catalog No:116906(ZBP1 antibody) at dilution of 1:50 (under 40x lens)