

PRMT3 antibody



Catalog Number: 114206

Product name

PRMT3 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

PRMT3 Rabbit Polyclonal antibody. Positive WB detected in human liver tissue. Positive IHC detected in human brain tissue. Observed molecular weight by Western-blot: 68 kDa

Preparation

This antibody was obtained by immunization of PRMT3 recombinant protein (Accession Number: NM_005788). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, IHC, WB


Dilutions

Recommended Dilution:

WB: 1:200-1:2000

IHC: 1:20-1:200


Validations



human liver tissue were subjected to SDS PAGE followed by western blot with Catalog No:114206 (PRMT3 antibody) at dilution of 1:600



Immunohistochemical of paraffin-embedded human brain using Catalog No:114206 (PRMT3 antibody) at dilution of 1:100 (under 10x lens)



Immunohistochemical of paraffin-embedded human brain using Catalog No:114206 (PRMT3 antibody) at dilution of 1:100 (under 40x lens)