

PAX3 antibody

Catalog Number: 113606

Product name

PAX3 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

PAX3 Rabbit Polyclonal antibody. Positive IP detected in mouse skin tissue. Positive WB detected in mouse skin tissue, HEK-293 cells. Observed molecular weight by Western-blot: 56-60 kDa

Preparation

This antibody was obtained by immunization of PAX3 recombinant protein (Accession Number: NM_181458). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IP


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IP: 1:200-1:2000

Validations



mouse skin tissue were subjected to SDS PAGE followed by western blot with Catalog No:113606(PAX3 antibody) at dilution of 1:500



IP Result of anti-PAX3 (IP:Catalog No:113606, 4ug; Detection:Catalog No:113606 1:500) with mouse skin tissue lysate 3200ug.