

RBM24 antibody

Catalog Number: 114606

Product name

RBM24 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

RBM24 Rabbit Polyclonal antibody. Positive IHC detected in human skeletal muscle tissue. Positive WB detected in HeLa cells.

Preparation

This antibody was obtained by immunization of RBM24 recombinant protein (Accession Number: NM_153020). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, IHC, WB


Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IHC: 1:20-1:200

Validations



HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:114606(RBM24 Antibody) at dilution of 1:600



Immunohistochemistry of paraffin-embedded human skeletal muscle tissue slide using Catalog No:114606(RBM24 Antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemistry of paraffin-embedded human skeletal muscle tissue slide using Catalog No:114606(RBM24 Antibody) at dilution of 1:50 (under 40x lens)