

SSBP1 antibody



Catalog Number: 115606

Product name

SSBP1 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

SSBP1 Rabbit Polyclonal antibody. Positive WB detected in A375 cells, mouse pancreas tissue, rat pancreas tissue. Observed molecular weight by Western-blot: 17 kDa

Preparation

This antibody was obtained by immunization of SSBP1 recombinant protein (Accession Number: NM_001256512). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:1000-1:10000

Validations



A375 cells were subjected to SDS PAGE followed by western blot with Catalog No:115606(SSBP1 antibody) at dilution of 1:500