

Human GALNT3 (His tag) recombinant protein



Catalog Number: 519606

General Information

Gene Name Synonym

Polypeptide N-acetylgalactosaminyltransferase 3, Polypeptide GalNAc transferase 3, GalNAc-T3, pp-GaNTase 3, Protein-UDP acetylgalactosaminyltransferase 3, UDP-GalNAc:polypeptide N-acetylgalactosaminyltransferase 3, HFTC, HHS

Protein Construction

A DNA sequence encoding the Human GALNT3 Gln38-Asp633 is expressed with a 6His tag at the C-terminus.

Organism

Human

Expression Host

Human Cells

QC Testing

Purity

>95% as determined by SDS-PAGE.

Endotoxin

<1.0 EU per μ g protein as determined by the LAL method.

Stability

Store for up to 6 months at -20°C

Formulation


Supplied as a $0.2\ \mu\text{m}$ filtered solution of 20mM PB, 150mM NaCl, pH7.4.

Usage Guide

Storage

Store at $<-20^{\circ}\text{C}$, minimize freeze-thaw cycles.

SDS-PAGE



Human GALNT3 (His tag) recombinant protein