

Anti-CD46 antibody



Catalog Number: 104306

Product name

Anti-CD46 antibody

Immunogen

[Human CD46 \(His Tag\) recombinant protein](#)

Specificity

Human CD46

Antibody description

Mouse monoclonal to CD46

Preparation

This antibody was produced from a hybridoma resulting from the fusion of a mouse myeloma with B cells obtained from a mouse immunized with purified, recombinant Human CD46 (rh CD46; EAW93471.1; Met 1-Asp 328) and conjugated with FITC under optimum conditions, the unreacted FITC was removed.

Formulation

Aqueous solution containing 0.5% BSA and 0.09% sodium azide

Storage

This antibody is stable for 12 months from date of receipt when stored at 2°C-8°C. Protected from prolonged exposure to light. Do not freeze !

Sodium azide is toxic to cells and should be disposed of properly. Flush with large volumes of water during disposal.

Clonality

Monoclonal

Ig Type

Mouse IgG1


Applications

FCM

Dilutions

FCM

Validations



CD46 Antibody (FITC), Mouse MAb, Flow cytometric

Profile of anti-CD46 reactivity on HeLa cells analyzed by flow cytometry.